

What causes changes to the Extracellular Matrix (ECM)?

With aging or during an injury, there are changes in the ECM both qualitatively and quantitatively¹. This affects the healing and recovery process from tendinopathy. ECMs are generally remodelled during an individual's life, but for tendons/ligaments this occurs at a slow rate. Parameters such as slow turnover, long-term post-translational modifications and extensive cell–matrix interactions are three aspects of ECM biology. Overuse and/or aging of tendon/ligament are major problems as this affects the mobility and with often slow rate of repair^{1-3,4,5}.


The oxygen consumption of tendons and ligaments is 7.5 times lower than that of skeletal muscles. The low metabolic rate and well-developed anaerobic energy-generation capacity are essential to carry loads and maintain tension for long periods, reducing the risk of ischemia and subsequent necrosis. However, a low metabolic rate results in slow healing after injury⁴.

Injecting into a tendon – penetrating a tendon sheath

What is a tendon sheath?

The entire tendon is covered by the epitenon, a fine, loose connective tissue sheath containing the vascular, lymphatic and nerve supply. More superficially, the epitenon is surrounded by paratenon, a loose areolar connective tissue consisting essentially of type I and type III collagen fibrils, some elastic fibrils, and an inner lining of synovial cells. Together, the paratenon and epitenon are sometimes called the peritendon³.

The epitenon is a white fibrous sheath surrounding a tendon. The sheath helps to reduce friction and provides structural support to the tendon. The fibrous consistency of the epitenon might create a form of resistance for an injection into the sheath. The epitenon contains a few fibroblast-like cells, blood and lymphatic vessels and a nerve supply. That is why injecting **into** a tendon sheath might result in pain.




Risk of rupture

There is a certain risk of causing further damage to the tendon when injecting into its sheath as this could result in tendon rupture⁶⁻⁸. The mechanics of such rupture could be associated with the sudden increase in volume in that small area, technique of injection, etc.

Why would an injection outside tendon sheath work?

The tendon and its sheath are surrounded by ECM which provides nutrients as well as support³. The tendon cells are embedded in a great amount of ECM which is secreted by the cells. It consists of protein fibres embedded in an amorphous mixture of huge protein-polysaccharide ("**proteoglycan**") molecules^{9,10}.

The contents of the ECM are dynamic and are able to flow freely between the tendon and the tendon sheath^{1,10}. Therefore, a peri-tendinous injection along the tendon sheath will be effective without causing the pain of penetrating the sheath¹¹. "The 'ground substance' of extracellular matrix is an amorphous gelatinous material. It is transparent, colourless, and fills the spaces between fibres and cells. It actually consists of large molecules called glycosaminoglycans (GAGs) which link together to form even larger molecules called proteoglycans. These molecules are very good at absorbing water, rather like a sponge, such that 90% of the extracellular matrix is made up of water. This means that the ECM is very good at resisting compressive forces."^{9,12}.




How does healing occur?

When a tendon is injured, it will naturally try to heal albeit slowly. This is known as the proliferative stage.

Proliferative Stage


The continued recruitment of fibroblasts and their rapid proliferation at the wound site are responsible for the synthesis of collagens, **proteoglycans**, and other components of the ECM. These components are initially arranged in a random manner within the ECM, which at this point is composed largely of type III collagen^{2,9}. An extensive blood vessel network is present, and the wound has a scar-like appearance. At the end of the proliferative stage, the repair tissue is highly cellular and contains relatively large amounts of water and an abundance of ECM components². Collagen type-III is gradually replaced by type-I collagen as the scar tissue matures^{1,11}.

Collagens and **proteoglycans** are components of the ECM critical to tendinogenesis^{1,12}.



Proteoglycans and Hyaluronic Acid

Proteoglycans are proteins that are heavily glycosylated. The basic proteoglycan unit consists of a "core protein" with one or more covalently attached glycosaminoglycan (GAG) chain(s). The backbone of a proteoglycan aggregate is a naturally occurring polymer called **Hyaluronic Acid (HA)**^{1,12}.



HA can also be found in the ECM. HA plays several important organizational roles in the ECM by binding with cells and other components through specific and nonspecific interactions. HA-binding proteins are constituents of the extracellular matrix, and stabilise its integrity¹³.

HA plays a fundamental role in the healing of an injured tendon.

How does STABHA work?

STABHA is a specifically manufactured soft tissue adapted biocompatible Hyaluronic Acid for soft tissue injuries (i.e. tendons and ligaments). STABHA is patented and is the only HA available for both the use in ligaments and tendons. Due to its superior purity and biocompatibility, STABHA is able to work efficiently without cause to worry for adverse effects that might be seen with other HA products.

STABHA is the supplementation of biocompatible HA to the injury site. This supplementation

ensures a "localised" delivery (via injection) to the injury site where HA is necessary in aiding the healing process.

Why more injections might be helpful?

The total normal turnover of HA in a human is 10-100mg/day¹³. Therefore, with one injection, it can be assumed that within 24 hours STABHA will be naturally eliminated by the body. Hence it is recommended for all three indications that there is a second injection of STABHA to ensure maintenance of biocompatible HA content at the injury site. It could be considered at more injections might be beneficial due to the rate of turnover during the recovery phase of the injury especially in older patients or patients with chronic tendinopathy.

References:

1. The pathogenesis of tendinopathy. A molecular perspective. Riley et al. *Rheumatology* 43 (2): 131–142.
2. Tendon: Biology, Biomechanics, Repair, Growth Factors, and Evolving Treatment Options. James et al. Review article, 2008 ASSH, published by Elsevier.
3. Tendon Injury and Tendinopathy: Healing and Repair. Sharma et al. *J. Bone Joint Surg. Am.* 87:187-202,2005.
4. Biology of tendon injury: healing, modeling and remodelling. Sharma et al. *J Musculoskeletal Neuronal Interact* 2006; 6(2):181-190.
5. The painful nonruptured tendon: clinical aspects. Khan et al. *Clin Sports Med* 22 (2003) 711–725.
6. Connecticut Center for Orthopedic Surgery. James T. Mazzara MD. www.orthoontheweb.com
7. Humpal Physical Therapy and Sports Medicine Centres. <http://www.humpalphysicaltherapy.com/Injuries-Conditions/Foot/Foot-Issues/Posterior-Tibial-Tendon-Problems-Patient-Guide/a~4995/article.html>
8. Common extensor tendon rupture following corticosteroid injection for lateral tendinosis of the elbow. Smith et al. *Br J Sports Med* 1999;33:423–425.
9. Rehabilitation of the Hand and Upper Extremity. Skirven et al. 6th edition published by Elsevier.
10. Regenerative Medicine and Biomaterials for the Repair of Connective Tissues. Archer et al. Published Jan 12, 2014 by Elsevier.
11. Intrinsic Healing Capacity and Tearing Process of Torn Supraspinatus Tendons: *In Situ* Hybridization Study of $\alpha 1(I)$ Procollagen mRNA. Hamada et al. *J. Orthop. Res.* Vol 15; No. 1: 1997
12. Histology Guide. Faculty of Biological Sciences, University of Leeds, United Kingdom
13. Hyaluronic acid (hyaluronan): a review. Necas et al. *Veterinarni Medicina*, 53, 2008 (8): 397–411.